



Prescribed for Progress™

Intelli-Fit Rotary Touchscreen User Manual



Corporate Headquarters

10601 W Belmont Ave, Franklin Park, IL 60131 • U.S.A.
847.288.3300 • FAX: 847.288.3703
Service phone number: 800.351.3737 (toll-free within U.S.A., Canada)
Global Website: www.lifefitness.com

International Offices

AMERICAS

North America

SCIFIT

5151 S. 110th E. Avenue
Tulsa, OK, 74146 U.S.A
Toll Free: +1 (800) 745 1373
Email: Service@scifit.com
Service Website: www.scifit.com/service-information

Brazil

Life Fitness Brasil

Av. Rebouças, 2315
Pinheiros
São Paulo, SP 05401-300
BRAZIL
SAC: 0800 773 8282 option 2
Telephone: +55 (11) 3095 5200 option 2
Service Email: suportebr@lifefitness.com
Sales/Marketing Email:
vendasbr@lifefitness.com

Latin America and Caribbean*

Life Fitness, LLC

10601 W Belmont Ave
Franklin Park, IL 60131 U.S.A.
Telephone: (847) 288 3300
Service Email: customersupport@lifefitness.com
Sales/Marketing Email:
commercialsales@lifefitness.com

EUROPE, MIDDLE EAST, and AFRICA (EMEA)

Netherlands and Luxemburg

Life Fitness Atlantic BV

Bijdorpplein 25-31
2992 LB Barendrecht
THE NETHERLANDS
Telephone: (+31) 180 646 666
Service Email: service.benelux@lifefitness.com
Sales/Marketing Email:
marketing.benelux@lifefitness.com

United Kingdom

Life Fitness UK LTD

Queen Adelaide
Ely, Cambs, CB7 4UB
Telephone: General Office (+44) 1353.666017
Customer Support (+44) 1353.665507
Service Email: uk.support@lifefitness.com
Sales/Marketing Email: life@lifefitness.com

Germany, Austria, and Switzerland

Life Fitness Europe GMBH

Neuhofweg 9
85716 Unterschleißheim
GERMANY
Telephone:
+49 (0) 89 / 31775166 Germany
+43 (0) 1 / 6157198 Austria
+41 (0) 848 / 000901 Switzerland
Service Email: kundendienst@lifefitness.com
Sales/Marketing Email: vertrieb@lifefitness.com

Spain

Life Fitness IBERIA

C/Frederic Mompou 5,1ª
08960 Sant Just Desvern Barcelona
SPAIN
Telephone: (+34) 93.672.4660
Service Email: servicio.tecnico@lifefitness.com
Sales/Marketing Email:
info.iberia@lifefitness.com

Belgium

Life Fitness Benelux NV

Parc Industriel de Petit-Rechain
4800 Verviers
BELGIUM
Telephone: (+32) 87 300 942
Service Email: service.benelux@lifefitness.com
Sales/Marketing Email:
marketing.benelux@lifefitness.com

All Other EMEA Countries and Distributor Business EMEA*

Life Fitness Atlantic BV

Bijdorpplein 25-31
2992 LB Barendrecht
THE NETHERLANDS
Telephone: (+31) 180 646 644
Service Email:
EMEAServiceSupport@lifefitness.com

ASIA PACIFIC (AP)

Japan

Life Fitness Japan, Ltd

4-17-33 Minami Aoyama 1F/B1F
Minato-ku - Tokyo 107-0062
Japan
Telephone: (+81) 0120.114.482
Fax: (+81) 03-5770-5059
Service Email: service.lfj@lifefitness.com
Sales/Marketing Email:
sales@lifefitnessjapan.com

Hong Kong

Life Fitness Asia Pacific LTD

26/F, Global Trade Square
21 Wong Chuk Hang Road
Wong Chuk Hang
Hong Kong
Telephone: (+852) 25756262
Fax: (+852) 25756894
Service Email: service.hk@lifefitness.com
Sales/Marketing Email:
hongkong.sales@lifefitness.com

All Other Asia Pacific countries and distributor business Asia Pacific*

Life Fitness Asia Pacific LTD

26/F, Global Trade Square
21 Wong Chuk Hang Road
Wong Chuk Hang
Hong Kong
Telephone: (+852) 25756262
Fax: (+852) 25756894
Service Email: service.hk@lifefitness.com
Sales/Marketing Email:
hongkong.sales@lifefitness.com

*Also check www.SCIFIT.com for local representation or distributor/dealer

* Order online 24/7 at www.SCIFIT.com

User and Service Documents Link

<https://lifecyclefitness9512.zendesk.com/hc/en-us>

<https://lifecycleworld.lifecyclefitness.com>

Additional information is available online using the links above.

تتوفر معلومات إضافية على الإنترنت باستخدام الرابط أعلاه.

点击上面的链接可在线获取更多信息。

Flere oplysninger er tilgængelige online gennem linket ovenfor.

Bijkomende informatie is online beschikbaar via bovenstaande link.

Vous trouverez plus d'informations en ligne à l'aide du lien ci-dessus.

Zusätzliche Informationen finden Sie online über den oben angegebenen Link.

Ulteriori informazioni sono disponibili online utilizzando il link sopra riportato.

追加情報は上記リンクを使用してオンラインで利用可能です。

상기 링크를 통해 온라인에서 추가 정보를 볼 수 있습니다.

Informações adicionais estão disponíveis on-line, através do link acima.

Дополнительная информация доступна в интернете по ссылке, указанной выше.

Mediante el enlace anterior podrá acceder a información adicional en línea.

Ytterligare information finns online genom att använda länken ovan.

İnternet üzerinden daha fazla bilgi edinmek için yukarıdaki bağlantıyı kullanabilirsiniz.

هناك معلومات إضافية متاحة على الإنترنت باستخدام الرابط أعلاه.

Informazio osagarria eskuragarri dago goiko estekaren bidez.

Допълнителна информация можете да намерите онлайн, като използвате връзката по-горе.

Mitjançant l'enllaç anterior podreu accedir a informació addicional en línia.

使用上面的連結線上提供額外資訊。

Dodatne informacije možete pronaći na internetu sljedeći vezu iznad.

ከላይ የተቀመጠውን አገናኝ(ሊንክ) በመጠቀም መረጃዎችን አንላይን ያገኛሉ፡፡

Lisätietoja on saatavissa verkosta käyttämällä yllä olevaa linkkiä.

Wubetumi anya nsem afororo aka ho wo websait so denam asem a ewo atifi ho a wubemia so so.

Πρόσθετες πληροφορίες είναι διαθέσιμες online χρησιμοποιώντας το σύνδεσμο παραπάνω.

מידע נוסף הוסף לקבל באינטרנט באמצעות הקישור לעיל.

További információ elérhető online, a fenti hivatkozás segítségével.

Viðbótarupplýsingar eru fáanlegar á netinu með því að smella á tengilinn hér fyrir ofan.

Plus indicium per superum situm potes invenire.

മുകളിലുള്ള ലിങ്ക് ഉപയോഗിച്ച് ഓൺലൈനിൽ കൂടുതൽ വിവരങ്ങൾ ലഭ്യമാണ്.

Ytterligere informasjon er tilgjengelig på nettet via linken ovenfor.

Dodatkowe informacje są dostępne online pod powyższym odnośnikiem.

Informações adicionais estão disponíveis online a usar o link acima.

Informații suplimentare sunt disponibile online, utilizând link-ul de mai sus.

Dodatne informacije dostupne su na mreži putem gornjeg linka.

Ďalšie informácie sú dostupné online na vyššie uvedenom odkaze.

Table of Contents

Safety

Introduction.....	5
Safety Information.....	5

SERVICE AND TECHNICAL DATA

Maintenance.....	7
------------------	---

Console-Tactile Key Overview

Console/Tactile Key Overview.....	8
-----------------------------------	---

Machine Program Menus

Machine Program Menus.....	9
----------------------------	---

Data Entry Screen

Data Entry Screen.....	10
------------------------	----

Workout Programs

QUICK START.....	11
MANUAL.....	11
CONSTANT WORK.....	12
ISO-STRENGTH.....	12
HILLS.....	13
HEART RATE.....	13
FIT-QUIK.....	14
MORE.....	14
RANDOM.....	14
STRESS TEST.....	15
HEART FIT TEST.....	15
HEART FIT TRAINING.....	16
POWER FIT TEST.....	16
POWER FIT TRAINING.....	17

Additional Key Features

POP-UP MENU.....	18
CHANGE VIEW.....	18
INCREASE/DECREASE LEVEL.....	18
PAUSE.....	18
CONTINUE/STOP.....	19
FAN.....	19
POWER BURST.....	19

Set-Up Mode Access

Touchscreen Access.....	20
Tactile Access.....	20
Set-Up Mode.....	21
Rotary Set Up Mode Categories.....	22
Overlay.....	22
Model.....	22

Language.....	22
Units.....	22
Baud Rate.....	22
Serial Protocol.....	23
Watts Scalar.....	23
Watt Offset.....	23
Power Burst (Seconds).....	23
Utilization (Hours).....	23
LCD Brightness.....	23
LCD Dimming.....	23
Default Level.....	23
Maximum Workout Time.....	23
Fit-Quik Cardio Workout.....	23
Fit-Quik Watts Target (Watts).....	24
Fit-Quik Profile Cardio.....	24
Fit-Quik Time Cardio (Min-Sec).....	24
Fit-Quik Reverse.....	24
Fit-Quik RPM Target.....	24
Fit-Quik Challenge Time (Min-Sec).....	24
Fit-Quik Recovery (Min-Sec).....	24
Fit-Quik ISO Repeats.....	24
Reverse Banner.....	24
Reverse Time Banner.....	24
ISO Banner.....	25
ISO Fast Time.....	25
ISO Slow Time.....	25
Paused Timeout (Seconds).....	25
Program Timeout (Seconds).....	25

Heart Rate Transmitter

Heart Rate Transmitter.....	26
-----------------------------	----

Battery Check

Access.....	27
Verify.....	27
Exit.....	27

Information Screens

Main Information Screen.....	28
Product Information.....	28
Service Information.....	28
Intelli-Fit Information.....	29

Warranty

What is Covered.....	30
Who is Covered.....	30
How long is it covered.....	30
Who Pays Transportation and Insurance For Service.....	30
What We Will Do To Correct Covered Defects.....	30

What is not covered.....	30
What you must do.....	30
How to get replacement parts and service.....	31
Exclusive Warranty.....	31
Changes in Warranty Not Authorized.....	31
Effects of State Laws.....	31

Introduction

Thank you for your purchase of this SCIFIT exercise machine. We have incorporated the finest technology and ergonomic design into this machine to assist you in achieving your fitness goals.


This SCIFIT machine is designed for class “S” (commercial) and class “A” (high accuracy) and class “II a” (medical directive) applications. The intended uses are for general health and wellness, cardio & pulmonary rehabilitation, sports medicine, physical therapy, physical conditioning and strength training. SCIFIT equipment are suitable for placement in medical, wellness, fitness, and recreational facilities. The machine has 191 levels of console controlled resistance, which ensures that users of all fitness levels will be able to get a challenging but not overwhelming workout.


To review the latest guidelines and recommendations, visit the American College of Sports Medicine website: www.acsm.org


FCC Warning - Possible Radio / Television Interference


This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy, and if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.


Safety Information


 **DANGER:** To reduce the risk of electrical shock or injury from moving parts, always unplug the machine before cleaning or attempting any maintenance activity.


 **WARNING:** Health-related injuries may result from incorrect or excessive use of exercise equipment. It is **STRONGLY** recommended seeing a physician for a complete medical exam before undertaking an exercise program, particularly if the user has a family history of high blood pressure or heart disease, is over the age of 45, smokes, has high cholesterol, is obese, or has not exercised regularly in the past year. If, at any time while exercising, the user experiences faintness, dizziness, pain, or shortness of breath, he or she must stop well immediately.

 **WARNING:** To reduce the risk of burns, fire, electrical shock, or injury, it is imperative to connect each product to a properly grounded electrical outlet.

 **WARNING:** Heart rate monitoring systems may be inaccurate. Over exercising may result in serious injury or death. If you feel faint, stop exercising immediately.

 **WARNING:** Allow a distance of 40 in. (1m) between the widest part of the machine and other objects on either side. Provide 40 in (1m) between the front or rear of the SCIFIT machine to any other objects.

 **WARNING:** This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to <http://www.P65Warnings.ca.gov>

 **CAUTION:** Any changes or modifications to this equipment could void the product warranty.

- Life Fitness Family of Brands does not warrant nor guarantee that component parts used in the manufacture of products offered under the Life Fitness Family of Brands are latex-free. Users of these products must take all necessary precautions to prevent accidental contact that could lead to an adverse latex reaction.
- Read all warnings and instructions on each product prior to using your exercise machine. Save these instructions for future reference.
- If warnings are missing or damaged, please contact customer service immediately for replacement warning labels. Warning labels are shipped with every product and should be installed before the product is used. SCIFIT is not responsible for missing or damaged warning labels.
- Use the exercise machine only for its intended purposes as described in this manual. Do not use attachments not recommended by the manufacturer.

- Never insert objects into any opening in these products. If an object should drop inside, turn off the power, unplug the power cord from the outlet, and carefully retrieve it. If the item cannot be reached, contact customer support services..
- This machine is intended for indoor use only. Do not place in direct sunlight, high humidity or where water is present (i.e. pool or spa areas).
- This machine is not intended for use by persons with reduced physical sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the machine by a person responsible for their safety. Keep children under the age of 14 away from this machine.
- Never operate the machine with the air openings blocked. Keep air openings free of lint, hair, or any other obstructing material.
- Never place liquids of any type directly on the unit, except in an accessory tray or holder. Containers with lids are recommended. Do not use these products in bare feet. Always wear shoes. Wear shoes with rubber or high-traction soles. Do not use shoes with heels, leather soles, cleats or spikes. Make sure no stones are embedded in the soles.
- Keep all loose clothing, shoelaces, and towels away from moving parts.
- Do not reach into, or underneath, the unit or tip it on its side during operation.
- Do not allow other people to interfere in any way with the user or equipment during a workout.
- When the product is not in use, SCIFIT recommends unplugging the product. Disconnect from the electrical outlet when not in use, and before putting on or taking off parts.
- Free standing equipment shall be installed on a stable and leveled surface.
- Do not stand or sit on plastic shrouds.
- Use caution when mounting or dismounting the machine. Use the stationary handlebar whenever additional stability is required.

SERVICE AND TECHNICAL DATA


Maintenance

COMPONENT	ACTION	WHEN
Console	Clean with mild soap and water	Weekly

Console-Tactile Key Overview


Console/Tactile Key Overview

Console Overview



1	Water bottle holder
2	Active tactile key area
3	Fan vents
4	Personal item compartment
5	Fit-Key connecting port
6	Active touchscreen area

Tactile Key Overview





1	Selects workout programs or activates a special feature.
2	Use during the data entry process or during a workout to increase/decrease the workout level.
3	Begins a pre-programmed workout session with no data entry required.
4	Confirms the user's value selection during the data entry process.
5	Pauses the workout, then two choices appear: Continue or Stop.



Machine Program Menus

Machine Program Menus



Model Series: PRO1 * PRO2 * PRO1000 * SXT7000 * ISO Upright Bikes * ISO Recumbent Bikes * REX7000



Model Series: StepOne Recumbent Stepper



Model Series: Latitude Lateral Stability Trainer




Data Entry Screen

Data Entry Screen



1	Numerical keypad to directly enter data information.
2	Returns the screen to the previous category.
3	Removes the last digits entered from right to left when using the numerical keypad.
4	Alternate way to enter user data information.
5	Accepts the user information entered and advances to the next category.



1	Returns the screen to the previous category.
2	Accepts the user information entered and advances to the next category.
3	Adjust the user data information to be entered.





Workout Programs

QUICK START




Description: Allows the user to get on and go! The workout has no time limit and the intensity level is manually controlled.

Changing Intensity

Touch the  key or the  tactile key to activate the screen pop-up menu OR press the   tactile keys.

MANUAL



Description: Allows the user to manually set workout time, intensity level and user weight.


Screen 1: Enter the time length of workout.

Screen 2: Enter the level intensity (Max. 20).

Screen 3: Enter the user weight.

COMMENCE WORKOUT

CONSTANT WORK



Description: Workout resistance is automatically adjusted to compensate for any change in speed (RPM), thus the amount of work will remain constant at any given (Watts or METs) setting.

Screen 1: Select target to either show “Watts” or “METs”.


Screen 2: Enter the time length of workout.

Screen 3: Enter the target “Watts” or “METs” according to the selection on screen 1.

Screen 4: Enter the user weight.

COMMENCE WORKOUT

ISO-STRENGTH



Description: Allows the user to set a target RPM rate. Intensity increases if the user exceeds the target RPM.

Screen 1: Enter the time length of workout.

Screen 2: Enter the target RPM.

Screen 3: Enter the user weight.

COMMENCE WORKOUT

HILLS



Description: Allows the user to select from 7 profiles set at various pre-programmed intensities and intervals. Profiles available: Course, Sprints, Twin Peaks, Peak, Progressive, Multi-Peaks and Ramp.

Screen 1: Select 1 of 7 course profiles.

Screen 2: Enter the time length of workout.

Screen 3: Enter the level intensity (Max. 20).

Screen 4: Enter the user weight.

COMMENCE WORKOUT

HEART RATE



Description: Enables the user to set a target heart rate. The machine maintains the heart rate within a zone by altering the intensity during the workout.

ATTENTION: A heart rate transmitter chest strap is required. See page ?? for instructions. If you are not wearing a transmitter belt or the belt is not sending a signal correctly, after a short period of time a warning message will appear.

Screen 1: Enter the time length of workout.

NOTE: Minimum workout time is 12 min. including a 2 min. warm up.

Screen 2: Enter the target Heart Rate.

Screen 3: Enter the user weight.

COMMENCE WORKOUT



Description: Allows the user to select a two stage program of cardiovascular work with ISO-strength (Isokinetic). Work loads and time durations are preset within the Set-Up mode.

Screen 1: Enter the user weight.

COMMENCE WORKOUT

MORE



Description: Touching this key will show a second menu screen (machine dependent) of additional pre-programmed workouts.

RANDOM



Description: Random workout profiles each with its own level of difficulty are shown for the user to choose.

NOTE: The highest column represents the workout level selected. Columns below it are percentages of the level selected.

Screen 1: Select 1 of 7 course profiles.


Screen 2: Enter the time length of workout.

Screen 3: Enter the level intensity (Max. 20).

Screen 4: Enter the user weight.

COMMENCE WORKOUT

STRESS TEST



Description: A designed stress test to be used as a training program or used with other diagnostic devices (if applicable). The user sets a target wattage with timed intervals and preset wattage increments.

Screen 1: Enter the watts target.

Screen 2: Enter the the number of column increments.

Screen 3: Enter the the time duration for the total increments.

Screen 4: Enter the user weight.

COMMENCE WORKOUT

HEART FIT TEST



Description: This program is based on a pre-programmed workout designed to evaluate your heart rate recovery. The user can select from 1 to 4 levels.

ATTENTION: A heart rate transmitter chest strap is required. See page ?? for instructions. If you are not wearing a transmitter belt or the belt is not sending a signal correctly, after a short period of time a warning message will appear.

Screen 1: Select 1 of 4 level intensities (Beginner, Intermediate, Advanced and Sports Performance).

Screen 2: Enter the user weight.

COMMENCE WORKOUT

HEART FIT TRAINING



Description: Designed training programs based on the user's evaluated Heart Fit test. A training schedule with retesting is recommended. The user has three program profiles and four level intensities to choose from.

ATTENTION: A heart rate transmitter chest strap is required. See page ?? for instructions. If you are not wearing a transmitter belt or the belt is not sending a signal correctly, after a short period of time a warning message will appear.

Screen 1: Select 1 of 3 program profiles.

Screen 2: Select 1 of 4 level intensities (Beginner, Intermediate, Advanced and Sports Performance).

Screen 3: Enter the user weight.

COMMENCE WORKOUT


POWER FIT TEST



Description: This test is based on a pre-programmed workout designed to evaluate your "Power Drop Percentage", which gives a fitness indication of power and strength using watts as the main measurement.

Press the Power Fit Test icon on the screen.

COMMENCE WORKOUT



Description: Specifically designed preset watts training programs based off the user's Power Fit test evaluation. Three program profiles and four level intensities to choose from.

Screen 1: Select 1 of 3 program profiles.

Screen 2: Select 1 of 4 level intensities (Beginner, Intermediate, Advanced and Sports Performance).


Screen 3: Enter the user weight.

COMMENCE WORKOUT

Additional Key Features

POP-UP MENU

Once a workout commences a pop-up sub menu may be access to show additional keys by touching the **Menu** icon or pressing the



tactile key.



CHANGE VIEW



Description: Allows the user to change the viewing format of the workout profile and data displayed on the screen during the workout.

INCREASE/DECREASE LEVEL



Description: Allows the user to change the intensity level of the workout.

PAUSE



Description: Touching the [Pause] icon on the screen or pressing the [Pause/Clear] tactile key pauses the workout program at anytime. The workout display will then minimize and the user is shown the CONTINUE/STOP image at the bottom of the screen.

CONTINUE/STOP



CONTINUE: Restarts the workout program at the point where the user had paused it.

STOP: Touch the [Stop] icon or press the [Pause/Clear] tactile key to exit the workout and return to the main screen.

FAN



Activates two 3-speed cooling fans (Low, Medium, High) at anytime during the workout. Each touch of the icon will change the fan airflow speed.

POWER BURST




Description: Increases resistance for a short duration of time during the workout when activated. The time duration of the “Power Burst” can be increased or decreased in the “Set-up” mode.

Set-Up Mode Access

Touchscreen Access

Method 1: Touching the screen in sequence


1. Touch the upper left corner of the screen.
2. Then touch the left side of the “Quick Start” icon.
3. Finally touch the upper left corner of the screen again.



Tactile Access


Method 2: Touching the tactile keys in sequence

1. Press the [Enter] key.
2. Then press the [Two Dot] key.
3. Then press the [Five Dot] key.
4. Finally press the [Six Dot] key.




Set-Up Mode

1. Category Selection




2. **NOTE:** There will be times the screen will change to a data entry view to make adjustments. After making the adjustments touch the “Accept” key to return to the “Set-Up Mode” screen.

Changing Parameters



3. **NOTE:** [Accept] to save and return to the main menu screen. [Cancel] to not save and return to the main menu screen.

Saving The Changes



Rotary Set Up Mode Categories

When the “Set-Up Mode” is accessed it will show categories for both rotary and treadmill. The following information will only list categories concerning the rotary machines.

Overlay

Default: SET AT FACTORY

Indicates the type of key array is selected for the membrane switch attached to the console.

Model

Default: MACHINE DEPENDENT

Configures the type of machine the console has been mounted to.

Language

Default: ENGLISH

Allows the user to select the language of their choice to be displayed on the screen.

Units

Default: U.S.

Selects the one of two commonly used weight units: U.S. (Imperial) or Metric.

Baud Rate

Default: 9600

The console may be connected via the “COMM” port on the back of the console to an external monitoring device that can be used to track a user’s progress. The baud rate determines the speed of data transfer between the console and the device. There are six different speed options to select from: 4800, 9600, 19200, 38400, 57600 and 115200.

Serial Protocol

Default: CSAFE

This gives the user a choice between “CSAFE” or “ERGOLINE” protocol.

Watts Scalar

Default: 1.000

This is a multiplication factor to change the power (Watts) of the machine. The machine will be harder to operate if the ‘Watt Scalar’ is greater than 1 and easier to operate if it is less than 1

Watt Offset

Default: 0.000

For machine calibration use only. **(Not for customer use)**

Power Burst (Seconds)

Default: 00:08

This increases resistance for a short duration of time during the workout when activated. The duration length of the “Power Burst” can be increased or decreased here.

Utilization (Hours)

Default: 0.00

This tracks the amount of time (in hours) that the machine has been in use.

NOTE: Hours will not reset when a software upgrade is done.

LCD Brightness

Default: 14

Adjusts the brightness of the LCD back light. Range is 8-14

LCD Dimming

Default: Disabled

When enabled activates the LCD screen to dim when RPM’s are 25 or lower.

Default Level

Default: 1

Controls the starting resistance level when a workout begins. If the level number is changed in the set-up mode, it will appear in the “Level” category during the data entry mode of the Manual, Quick Start, Hills and Random programs.

Maximum Workout Time

Default: 99:00

Sets a workout time limit on programs Manual, Hills, ISO-Strength, Heart Rate and Constant Work. No workout time will be allowed to surpass the limit during the data entry mode.

Fit-Quik Cardio Workout

Default: Manual

Allows the Fit-Quik program to be preset to 1 of 3 workout programs: Manual, Hills or Constant Workout.

Fit-Quik Watts Target (Watts)

Default: 20

A preset wattage target is selected here when the Fit-Quik program is set for Constant Work.

Fit-Quik Profile Cardio

Default: Course

When the Fit-Quik program is preset for the Hills program the user may select 1 of 7 workout contours: Course, Sprint, Twin peaks, Peak, Progressive, Multi-Peaks and Ramp.

Fit-Quik Time Cardio (Min-Sec)

Default: 5:00

Sets the time length of the 1st stage workout for the program before initiating the 2nd stage.

Fit-Quik Reverse

Default: Disabled

Sets the option to display the phrase “Change Direction” on the screen during the 1st stage workout for those machines that are bi-directional compatible.

Fit-Quik RPM Target

Default: Machine Dependent

Sets the RPM target for the 2nd stage.

Fit-Quik Challenge Time (Min-Sec)

Default: 0:10

Allows the option to increase or decrease the time length of the ISO portion stage.

Fit-Quik Recovery (Min-Sec)

Default: 00:35

Allows the option to increase or decrease the time length of the Constant Work program portion.

Fit-Quik ISO Repeats

Defaults: 4

Allows the option to set how many times the ISO-Strength portion of the 2nd stage is repeated.

Reverse Banner

Default: Disabled

Allows the option to display the phrase “Change Direction” for bi-directional compatible machines.

Affects the following programs: Manual, Constant Work, Hills and Random.

Reverse Time Banner

Default: 00:30

Sets the time interval when the phrase “Change Direction” appears on the screen.

The phrase “Change Direction” will appear every 30 seconds as set by the default.

ISO Banner

Default: Disabled

Allows the option to enable or disable the banners “Fast Forward”, “Slow” and Reverse Fast” to be shown during an ISO-Strength workout.

ISO Fast Time

Default: 00:15

Sets the time interval before the phrases “Fast Forward” and “Reverse Fast” appear on the screen during an ISO-Strength workout.

Default: Every 15 seconds.

ISO Slow Time

Default: 00:15

Sets the time interval before the phrase “Slow” appears on the screen during an ISO-Strength workout. Default: Every 15 seconds

Paused Timeout (Seconds)

Default: 00:60

The amount of time the screen will stay lit after the [Pause] key is pressed, the user stops pedaling and the brake stops rotating. This time can be adjusted here.

NOTE: Note: Paused timeout range is 0 - 180 seconds.

Program Timeout (Seconds)

Default: 00:60

During data entry if no input is done within 60 seconds it will return to the main menu.

Heart Rate Transmitter

Heart Rate Transmitter




WARNING: Heart Rate monitoring systems may be inaccurate. If you feel faint, stop exercising immediately. Over exercising may result in serious injury or death.



CAUTION: Always follow manufacturers instructions for proper care and use. Failure to follow manufacturers instructions may result in injury to yourself.

ATTENTION: The wireless chest heart rate strap is required when selecting the Heart Rate, Heart Fit or Heart Fit training programs. The wireless strap will not work when worn over clothing. **Skin contact is necessary**



Calculating Target Heart Rate Zone

Your approximate Maximum Heart Rate (MHR) is equal to 220 minus your age. The upper limit of your target zone is equal to 0.85 times your MHR. The lower limit of your target zone is equal to 0.60 times your MHR.

Example

Your age is 40:

$$220 - 40 = 180$$

$$180 \times .85 = 153 \text{ (Upper target zone limit)}$$

$$180 \times .60 = 108 \text{ (Lower target zone limit)}$$


These values are based upon averages, always consult your physician to establish the proper heart rate zone for your individual health condition.

Battery Check

Access


Touching the screen in sequence

1. Touch the left side of the “Quick Start” icon.
2. Then touch the upper left corner of the screen.
3. Finally touch the left side of the “Quick Start” icon again.




Verify

If the battery voltage is below 10.5V, then replace the battery.



Exit


Press and hold the [1 Dot] tactile key, then, press the [Pause/Clear] tactile key and release both keys.



Information Screens

Main Information Screen


Touch the [More] key on the main menu screen until the information screen is shown (see below), then touch the information icon you wish to access.



Product Information

This screen will show the following:

1. Domestic website.
2. Touch the 'Safety Advice' icon to display additional recommended safety information.



Service Information

The screen will show the following:


1. Domestic website.
2. Current console software installed.



Intelli-Fit Information

This screen will show the following:

1. A description summary of the console programs available (Machine dependent).



Warranty

What is Covered

This Life Fitness Family of Brands commercial exercise equipment product is warranted to be free of all defects in material and workmanship.

Who is Covered

The original purchaser or any person receiving the product as a gift from the original purchaser.

How long is it covered

All coverage is provided by specific Product according to the guidelines listed on the chart below.

Item	6 Months	1 Year	2 Years	3 Years	5 Years
Frame					X
Electrical				X	
Mechanical				X	
Labor		X			
Cranks/Drive Syst. (Latitude Only)			X		
Upholstery, Seats, Grips, Walk belts, Decks		X			
Battery		X			
High Wear Items (Heartrate Straps, Pedal Straps)	X				

Who Pays Transportation and Insurance For Service

If the Product or any covered part must be returned to a service facility for repairs, We, Life Fitness Family of Brands, will pay all transportation and insurance charges for the first year. You are responsible for transportation and insurance charge after the first year.

What We Will Do To Correct Covered Defects

We will ship to you any new or rebuilt replacement part or component, or at our option, replace the Product. Such replacement parts are warranted for the remaining portion of the original warranty period.

What is not covered

Any failures or damage caused by unauthorized service, misuse, accident, negligence, improper assembly or installation, debris resulting from any construction activities in the Product's environment, rust or corrosion as a result of the Product's location, alterations or modifications without our written authorization or by failure on your part to use, operate and maintain the Product as set out in your Operation Manual ("Manual").

One type RJ45 interchangeable Network ready and Fitness Entertainment port is supplied with the Product. This port complies with the FitLinxx CSAFE specification dated August 4, 2004 stating: 4.75 VDC to 10 VDC; maximum current of 85 mA. Any Product damage caused by a load exceeding this FitLinxx CSAFE specification is not covered by warranty.

All terms of this warranty are void if this product is moved beyond the continental borders of the United States of America (excluding Alaska, Hawaii and Canada) and are then subject to the terms provided by that country's local authorized *SCIFIT* representative.

Detailed warranty guidelines are available at: <http://www.scifit.com/service-information/>.

What you must do

Register your Product online at <http://www.scifit.com/warranty-registration-2/>.

Retain proof of purchase; use, operate and maintain the Product as specified in the Manual; notify Customer Support Services of any defect within 10 days after discovery of the defect; if instructed, return any defective part for replacement or, if necessary, the entire Product for repair.

How to get replacement parts and service

Refer to the front of this manual for your local service contact information. Reference your name, address and the serial number of your Product. They will tell you how to get a replacement part, or, if necessary, arrange for service where your Product is located.

Exclusive Warranty

THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES OF ANY KIND EITHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, AND ALL OTHER OBLIGATIONS OR LIABILITIES ON OUR PART. We neither assume nor authorize any person to assure for us any other obligation or liability concerning the sale of this Product. Under no circumstances shall we be liable under this warranty, or otherwise, of any damage to any person or property, including any lost profits or lost savings, for any special, indirect, secondary, incidental or consequential damages of any nature arising out of the use of or inability to use this Product. Some states do not allow the exclusion or limitation of implied warranties or of liability for incidental or consequential damages, so the above limitations or exclusions may not apply to you. Warranty coverages and terms may differ outside the United States. Please contact the Life Fitness office servicing your country (contact information found at the front of this manual) or visit the applicable local Life Fitness website to receive the specific warranty information for your country.

Changes in Warranty Not Authorized

No one is authorized to change, modify or extend the terms of this limited warranty.

Effects of State Laws

This warranty gives you specific legal rights, and you may have other rights which vary from state to state and country by country.



Prescribed for Progress™

A6802C
2023

www.scifit.com